

# Reliable long-term storage solution

## CryoPlus sample storage systems

Thermo Scientific™ CryoPlus™ systems provide the perfect combination of liquid nitrogen storage reliability and microprocessor technology. Storage capacities of up to 39,000 2.0 mL vials (8,450 to 65,910 Thermo Scientific™ Nunc™ Cryobank™ vials) make the most of your laboratory storage space.



### Storage system features

- Top-mounted control panel allows easy access to the unit's microprocessor controller for programming
- All lids incorporate proprietary foamed-in-place, high-density urethane insulation with dedicated vent; two independent, flexible gaskets reduce the migration of moisture into the chamber and reduce icing
- 24 tricolor LEDs continuously display actual liquid nitrogen level and high-level/low-level setpoints
- Remote alarm contacts for in-house remote alarms or connection to a Thermo Scientific™ Sensaphone™ Telephone Dialing System
- Counterbalanced UL 61010 third edition–compliant lid with 100% clearance allows each entry into the storage area for retrieval of samples, while minimizing sample exposure time and maintaining a safe work environment
- Heavy-duty casters are provided standard
- Incorporates a temperature sleeve as a standard feature, which provides colder temperatures at the top and more efficient operation in vapor phase
- Storage tanks are vacuum-insulated, with stainless steel interior
- Exterior cabinet is constructed of 18-gauge, cold roll steel with powder-coat paint for durability and a high-quality uniform finish; salt spray tests exceed 1,000 hours per ASTM Standard 117B-85
- All models come standard with LN<sub>2</sub> transfer hose


### For maximum product security, CryoPlus systems incorporate three different solutions

- A convenient front-mounted cabinet key lock for sample security
- A keyed lock on the microprocessor for temperature and alarm setpoint security
- A rear-mounted, recessed power switch to prevent accidental power shut-off

### A complete sample storage offering for vapor or liquid phase

- Square and vertical racks for vial storage
- Frames and canisters for bag storage
- Specialty racks for cane storage
- Gas bypass assembly delivers warm nitrogen to atmosphere until the control probe reaches set temperature. The bypass valve closes and the fill valve is energized, preventing the storage area from being subjected to warm nitrogen gas. This is also required if the CryoPlus unit is connected to a manifold system to prevent LN<sub>2</sub> source alarms.

## CryoPlus system information panel



### CryoPlus sample storage systems

Model	Cat. No.	LN <sub>2</sub> capacity (liters)	Vial capacity 1.2–2.0 mL (CryoBank)	Electrical	Plug type	Inner tank diameter in. (cm)	Exterior dimensions H x W x D in. (cm)	Shipping weight	Regulatory listings	LN <sub>2</sub> and utility connection
CryoPlus 1 system	7400	90	6,318 <sup>1</sup> (8,450 <sup>2,3</sup> )	115v/60 Hz	NEMA 5-15P	16 (40.6)	41 x 21.5 x 26 (104.1 x 54.6 x 66)	260 lbs/ 117.9kg	CFDA, cULus and CE	All models require 22 psi low pressure supply tank and are supplied with (1) 0.5 in. 45° flare stainless steel flexible transfer hose
	7401			208-230v/50 Hz	CEE 7/7					
CryoPlus 2 system	7402	200	13,500 <sup>2</sup> (22,815 <sup>2,3</sup> )	115v/60 Hz	NEMA 5-15P	24 (61)	41 x 28.5 x 34 (104.1 x 72.4 x 86.4)	480 lbs/ 217.7 kg		
	7403			208-230v/50 Hz	CEE 7/7					
CryoPlus 3 system	7404	340	23,000 <sup>2</sup> (38,870 <sup>2,3</sup> )	115v/60 Hz	NEMA 5-15P	31 (78.7)	41 x 34.5 x 41.5 (104.1 x 87.6 x 105.4)	550 lbs/ 249 kg		
	7405			208-230v/50 Hz	CEE 7/7					
CryoPlus 4 system	7406	552	39,000 <sup>2</sup> (65,910 <sup>2,3</sup> )	115v/60 Hz	NEMA 5-15P	39.5 (110.3)	47 x 43.5 x 50 (119.4 x 110.5 x 127)	1070 lbs/ 484.3 kg		
	7407			208-230v/50 Hz	CEE 7/7					


1 Using arrowhead racks no. 4000044.

2 Using vertical racks no. 4000012 and dense storage tube boxes.

3 Using Thermo Scientific™ Nunc™ Cryobank™ dense storage boxes (196-cell) and 1 mL Cryobank tubes.

### Thermo Scientific™ CryoPlus™ vapor phase starter kits

Model	Cat. No.	LN <sub>2</sub> capacity (liters)	Electrical	Description
CryoPlus 1 Vapor Phase Starter Kit	7400RAKHGBR	90	120V 50/60Hz	Package includes CryoPlus 1 (model 7400) and a hot gas bypass option, 4 storage rack (rack 4000001), and a 3.5" riser, UL
	7401RAKHGBR		230V 50/60Hz	Package includes CryoPlus 1 (model 7401) and a hot gas bypass option, 4 storage rack (rack 4000001), and a 3.5" riser, CE
CryoPlus 2 Vapor Phase Starter Kit	7402RAKHGBR	200	120V 50/60Hz	Package includes CryoPlus 2 (model 7402) and a hot gas bypass option, 10 storage rack (rack 4000001), and a 3.5" riser, UL
	7403RAKHGBR		230V 50/60Hz	Package includes CryoPlus 2 (model 7403) and a hot gas bypass option, 10 storage rack (rack 4000001), and a 3.5" riser, CE
CryoPlus 3 Vapor Phase Starter Kit	7404RAKHGBR	340	120V 50/60Hz	Package includes CryoPlus 3 (model 7404) and a hot gas bypass option, 17 storage rack (rack 4000001), and a 3.5" riser, UL
	7405RAKHGBR		230V 50/60Hz	Package includes CryoPlus 3 (model 7405) and a hot gas bypass option, 17 storage rack (rack 4000001), and a 3.5" riser, CE
CryoPlus 4 Vapor Phase Starter Kit	7406RAKHGBR	552	120V 50/60Hz	Package includes CryoPlus 4 (model 7406) and a hot gas bypass option, 29 storage rack (rack 4000001), and a 3.5" riser, UL
	7407RAKHGBR		230V 50/60Hz	Package includes CryoPlus 4 (model 7407) and a hot gas bypass option, 29 storage rack (rack 4000001), and a 3.5" riser, CE




### CryoPlus rack selection guide

Cat. No. and description		Dimensions H x W x D in. (cm)	Storage	CryoPlus 1 system	CryoPlus 2 system	CryoPlus 3 system	CryoPlus 4 system
<b>Square racks for liquid phase storage</b>							
	<b>4000006*</b>	27.4 x 5.5 x 5.63 (69.6 x 14 x 14.3)	Racks per vessel	4	10	17	29
	Rack includes cardboard boxes and cardboard dividers for 1.5/2 mL vials		Boxes per rack	13	13	13	13
	Vials per box		100	100	100	100	
	Vials per vessel		5,200	13,000	22,100	37,700	
	<b>4000007</b>	27.4 x 5.5 x 5.63 (69.6 x 14 x 14.3)	Racks per vessel	4	10	17	29
	Rack includes stainless steel boxes and cardboard dividers for 1.5/2 mL vials		Boxes per rack	13	13	13	13
	Vials per box		100	100	100	100	
	Vials per vessel		5,200	13,000	22,100	37,700	
	<b>4000042</b>	26 x 5.5 x 5.63 (66 x 14 x 14.3)	Racks per vessel	4	10	17	29
	Rack includes cardboard boxes and cardboard dividers for 1.5/2 mL vials		Boxes per rack	12	12	12	12
	Vials per box		100	100	100	100	
	Vials per vessel		4,800	12,000	20,400	34,800	
	<b>4000008</b>	25.4 x 5.5 x 5.63 (64.5 x 14 x 14.3)	Racks per vessel	4	10	17	29
	Rack includes cardboard boxes and cardboard dividers for 4 mL vials		Boxes per rack	8	8	8	8
	Vials per box		100	100	100	100	
	Vials per vessel		3,200	8,000	13,600	23,200	
	<b>4000009</b>	25.4 x 5.5 x 5.63 (64.5 x 14 x 14.3)	Racks per vessel	4	10	17	29
	Rack includes stainless steel boxes and cardboard dividers for 4 mL vials		Boxes per rack	8	8	8	8
	Vials per box		100	100	100	100	
	Vials per vessel		3,200	8,000	13,600	23,200	
	<b>4000010</b>	26.6 x 5.5 x 5.63 (67.6 x 14 x 14.3)	Racks per vessel	4	10	17	29
	Rack includes cardboard boxes and cardboard dividers for 5 mL vials		Boxes per rack	7	7	7	7
	Vials per box		100	100	100	100	
	Vials per vessel		2,800	7,000	11,900	20,300	

\* Rack no. 4000006 will not accommodate polycarbonate boxes; use rack no. 4000042 instead.


**CryoPlus rack selection**

Cat. No. and description		Dimensions H x W x D in. (cm)	Storage	CryoPlus 1 system	CryoPlus 2 system	CryoPlus 3 system	CryoPlus 4 system
<b>Square racks for vapor phase storage</b>							
	<b>4000001</b>	23.9 x 5.5 x 5.63 (60.7 x 14 x 14.3)	Racks per vessel	4	10	17	29
	Rack includes cardboard boxes and cardboard dividers for 1.5/2 mL vials		Boxes per rack	11	11	11	11
			Vials per box	100	100	100	100
			Vials per vessel	4,400	11,000	18,700	31,900
	<b>4000002</b>	23.9 x 5.5 x 5.63 (60.7 x 14 x 14.3)	Racks per vessel	4	10	17	29
	Rack includes stainless steel boxes for 1.5/2 mL vials		Boxes per rack	11	11	11	11
			Vials per box	100	100	100	100
			Vials per vessel	4,400	11,000	18,700	31,900
	<b>4000003</b>	22.3 x 5.5 x 5.63 (56.6 x 14 x 14.3)	Racks per vessel	4	10	17	29
	Rack includes cardboard boxes and cardboard dividers for 4 mL vials		Boxes per rack	7	7	7	7
			Vials per box	100	100	100	100
			Vials per vessel	2,800	7,000	11,900	20,300
	<b>4000004</b>	22.3 x 5.5 x 5.63 (56.6 x 14 x 14.3)	Racks per vessel	4	10	17	29
	Rack includes stainless steel boxes and cardboard dividers for 4 mL vials		Boxes per rack	7	7	7	7
			Vials per box	100	100	100	100
			Vials per vessel	2,800	7,000	11,900	20,300
	<b>4000005</b>	22.8 x 5.5 x 5.63 (57.9 x 14 x 14.3)	Racks per vessel	4	10	17	29
	Rack includes cardboard boxes and cardboard dividers for 5 mL vials		Boxes per rack	6	6	6	6
			Vials per box	100	100	100	100
			Vials per vessel	2,400	6,000	10,200	17,400
<b>Specialty racks</b>							
	<b>4000379</b>	23.9 x 5.5 x 5.63 (60.7 x 14 x 14.3)	Racks per vessel	4	10	17	29
	SUC-1 rack for cane storage (includes canisters #4000176)		Canisters per rack	8	8	8	8
			Canes per canister	25	25	25	25
			Canes per vessel	800	2,000	3,400	5,800
<b>Platform risers required for racks in vapor phase storage</b>							
	3.5" riser for use with 4000001, 4000002, and 4000379 racks			<b>4000060</b>	<b>4000060</b>	<b>4000063</b>	<b>4000066</b>
	4" riser for use with 4000005 racks			<b>4000056</b>	<b>4000062</b>	<b>4000065</b>	<b>4000068</b>
	5" riser for use with 4000003 and 4000004 racks			<b>4000061</b>	<b>4000061</b>	<b>4000064</b>	<b>4000067</b>



### Arrowhead racks for liquid phase storage

Cat. No. and description	Dimensions H x W x D in. (cm)	Storage	CryoPlus 1 system
 <p data-bbox="266 352 971 537">4000044 Rack includes cardboard boxes and cardboard dividers for 1.5/2 mL vials</p>	28.1 x 7.8 x 7.5 (71.4 x 19.8 x 19.1)	Racks per vessel	6
		Boxes per rack	13
		Vials per box	81
		Vials per vessel	6,318
 <p data-bbox="266 552 971 737">4000112 Rack includes cardboard boxes and cardboard dividers for 1.5/2 mL vials</p>	26 x 7.8 x 7.5 (66 x 19.8 x 19.1)	Racks per vessel	6
		Boxes per rack	12
		Vials per box	81
		Vials per vessel	5,832
 <p data-bbox="266 751 971 936">4000113 Rack includes stainless steel boxes and cardboard dividers for 1.5/2 mL vials</p>	26 x 7.8 x 7.5 (66 x 19.8 x 19.1)	Racks per vessel	6
		Boxes per rack	12
		Vials per box	81
		Vials per vessel	5,832
 <p data-bbox="266 942 971 1136">4000104 Rack includes cardboard boxes and cardboard dividers for 4 mL vials</p>	25.4 x 7.8 x 7.5 (64.5x 19.8 x 19.1)	Racks per vessel	6
		Boxes per rack	8
		Vials per box	81
		Vials per vessel	3,888
 <p data-bbox="266 1142 971 1335">4000105 Rack includes stainless steel boxes and cardboard dividers for 4 mL vials</p>	25.4 x 7.8 x 7.5 (64.5x 19.8 x 19.1)	Racks per vessel	6
		Boxes per rack	8
		Vials per box	81
		Vials per vessel	3,888
 <p data-bbox="266 1341 971 1535">4000043 Rack includes cardboard boxes and cardboard dividers for 5 mL vials</p>	26.5 x 7.8 x 7.5 (67.3 x 19.8 x 19)	Racks per vessel	6
		Boxes per rack	7
		Vials per box	81
		Vials per vessel	3,402



### Arrowhead racks for vapor phase storage


Cat. No. and description		Dimensions H x W x D in. (cm)	Storage	CryoPlus 1 system
	4000110	23.8 x 7.8 x 7.5 (60.5 x 19.8 x 19.1)	Racks per vessel	6
	Rack includes cardboard boxes and cardboard dividers for 1.5/2 mL vials		Boxes per rack	11
			Vials per box	81
			Vials per vessel	5,346
	4000111	23.8 x 7.8 x 7.5 (60.5 x 19.8 x 19.1)	Racks per vessel	6
	Rack includes stainless steel boxes and cardboard dividers for 1.5/2 mL vials		Boxes per rack	11
			Vials per box	81
			Vials per vessel	5,346
	4000108	21.6 x 7.8 x 7.5 (54.9 x 19.8 x 19.1)	Racks per vessel	6
	Rack includes cardboard boxes and cardboard dividers for 2 mL vials		Boxes per rack	10
			Vials per box	81
			Vials per vessel	4,860

### Platform risers required for racks in vapor phase storage

	Description of riser and racks for intended use	CryoPlus 1 system
	3.5 in. riser for use with 4000110 and 4000111 racks	4000060
	5 in. riser for use with 4000108 racks	4000061

### Arrowhead boxes and dividers for vapor phase storage

	Cat. No.	Description
	<b>Boxes</b>	
	4000193	Arrowhead cardboard box (with drain hole) includes 81-cell divider for use with 1.5/2 mL vials
	4000194	Arrowhead stainless steel box (with drain hole) does not include 81-cell divider for use with 1.5/2 mL vials
	4000195	Arrowhead cardboard box (with drain hole) includes 81-cell divider for use with 4/5 mL vials
	4000196	Arrowhead stainless steel box (with drain hole) does not include 81-cell divider for use with 4 mL vials
	<b>Cell dividers</b>	
4000181	81-cell divider for use with all arrowhead boxes	



### Vertical racks for liquid phase storage

Cat. No. and description	Dimensions H x W x D in. (cm)	Storage	CryoPlus 1 system	CryoPlus 2 system	CryoPlus 3 system	CryoPlus 4 system	
<b>Vertical rack for liquid phase storage</b>							
	4000012	27.6 x 2.3 x 5.7 (70.1 x 5.8 x 14.5)	Racks per vessel	10	27	46	78
	Rack includes cardboard boxes and cardboard dividers for 1.5/2 mL vials		Boxes per rack	5	5	5	5
	Vials per box		100	100	100	100	
	Vials per vessel		5,000	13,500	23,000	39,000	

<b>Vertical rack for vapor phase storage</b>							
	4000011	22.2 x 2.3 x 5.7 (56.4 x 5.8 x 14.5)	Racks per vessel	10	27	46	78
	Rack includes cardboard boxes and cardboard dividers for 1.5/2 mL vials		Boxes per rack	4	4	4	4
	Vials per box		100	100	100	100	
	Vials per vessel		4,000	10,800	18,400	31,200	

<b>Platform risers required for racks in vapor phase storage</b>							
<b>Description of riser and racks for intended use</b>							
	5 in. riser for use with 4000011 racks			4000061	4000061	4000064	4000067

### Boxes and dividers

Cat. No.	Description
 4000014	Cardboard box with 100-cell divider for use with 1.5/2 mL vials (does have drain hole)
820010	Polycarbonate box with 100-cell divider for use with 1.5/2 mL vials with internal threaded caps (does have drain hole)
820011	Polycarbonate box with 81-cell divider for use with 1.5/2 mL vials (does have drain hole)

### Lid locking pins

Cat. No.	Description
1890632	Factory-installed lid-locking pins (allows customers to securely lock the lid in the open position as a secondary user-safety feature)*
1890633	Certified technician (field-installed) lid-locking pins (allows customers to securely lock the lid in the open position as a secondary user-safety feature)

\* Lid-locking pins are not available as factory-installed items when a CryoPlus vapor-phase starter SKU is ordered.

### Frames and canisters for liquid phase storage

Cat. No. and description		Dimensions H x W x D in. (cm)	Storage	CryoPlus 1 system	CryoPlus 2 system	CryoPlus 3 system	CryoPlus 4 system
	4000670	27.4 x 0.6 x 4.6 (69.6 x 1.5 x 11.7)	Frames per vessel	28	72	124	204
	Frame for 4000610 canister		Canisters per frame	7	7	7	7
			Canisters per vessel	196	504	868	1,428
	1900585	27.4 x 0.6 x 3.8 (69.6 x 1.5 x 9.7)	Frames per vessel	34	100	184	314
	Frame for 1950831 canister		Canisters per frame	7	7	7	7
			Canisters per vessel	238	700	1,288	2,198
	4000371	23.7 x 0.6 x 6.5 (60.2 x 1.5 x 16.5)	Frames per vessel	28	68	108	186
	Frame for 4000335 canister		Canisters per frame	4	4	4	4
			Canisters per vessel	112	272	432	744
	4000372	23.8 x 0.7 x 8.2 (60.5 x 1.8 x 20.8)	Frames per vessel	22	50	90	154
	Frame for 4000336 canister		Canisters per frame	4	4	4	4
			Canisters per vessel	88	200	360	616
	4000368	26.4 x 0.9 x 6.1 (67.1 x 2.3 x 15.5)	Frames per vessel	19	46	72	126
	Frame for 4000332 and 4000356 canisters		Canisters per frame	4	4	4	4
			Canisters per vessel	76	184	288	504
	4000369	26.4 x 1 x 10 (67.1 x 2.5 x 25.4)	Frames per vessel	10	28	48	84
	Frame for 4000333 and 4000357 canisters		Canisters per frame	4	4	4	4
			Canisters per vessel	40	112	192	336

### Frames and canisters for vapor phase storage

Cat. No. and description		Dimensions H x W x D in. (cm)	Storage	CryoPlus 1 system	CryoPlus 2 system	CryoPlus 3 system	CryoPlus 4 system
	4000660	23.5 x 0.6 x 4.6 (59.7 x 1.5 x 11.7)	Frames per vessel	28	72	124	204
	Frame for 4000610 canister		Canisters per frame	6	6	6	6
			Canisters per vessel	168	432	744	1,224
	1900584	23.5 x 0.6 x 3.8 (59.7 x 1.5 x 9.7)	Frames per vessel	34	100	184	314
	Frame for 1950831 canister		Canisters per frame	6	6	6	6
			Canisters per vessel	204	600	1,104	1,884
	4000371	23.7 x 0.6 x 6.5 (60.2 x 1.5 x 16.5)	Frames per vessel	28	68	108	186
	Frame for 4000335 canister		Canisters per frame	4	4	4	4
			Canisters per vessel	112	272	432	744
	4000372	23.8 x 0.7 x 8.2 (60.5 x 1.8 x 20.8)	Frames per vessel	22	50	90	154
	Frame for 4000336 canister		Canisters per frame	4	4	4	4
			Canisters per vessel	88	200	360	616

## Platform risers required for frames in vapor phase storage


	Description of riser and frames for intended use	CryoPlus 1 system	CryoPlus 2 system	CryoPlus 3 system	CryoPlus 4 system
	3.5 in. riser for use with 4000371 and 4000372 frames	4000060	4000060	4000063	4000066
	4 in. riser for use with 4000660 frames	4000056	4000062	4000065	4000068

## Platform dividers

	Description of platform divider and frames for intended use	CryoPlus 1 system	CryoPlus 2 system	CryoPlus 3 system	CryoPlus 4 system
	Platform divider for use with 4000660 and 4000670 25 mL frames: 1900584 and 1900585	4000616	4000624	4000631	4000639
	Platform divider for use with 4000371 frames	N/A	4000075	4000050	4000053
	Platform divider for use with 4000372 frames	N/A	4000076	4000051	4000054
	Platform divider for use with 4000368 frames	N/A	N/A	4000047	4000052
	Platform divider for use with 4000369 frames	N/A	N/A	4000048	N/A

## Canisters

	Cat. No. and description	Dimensions H x W x D in. (cm)
	<b>1950831</b> Swing arm canister for 25 mL bag (OriGen CS 25 or equivalent)	3.6 x 3.915 x .375 (9.14 x 9.94 x 0.95)
	<b>4000610</b> Swing arm canister for 50 mL bag (Fenwal 4R9951 or OriGen CS 50)	3.7 x 6.3 x 0.5 (9.4 x 16 x 1.3)
	<b>4000335</b> Sliding canister for 250 mL bag (Fenwal 4R5461 or OriGen CS 250, CryoMACS 50 and 250)	5.5 x 7.6 x 0.4 (14 x 19.3 x 1)
	<b>4000336</b> Sliding canister for 500 mL bag (Fenwal 4R5462 or OriGen CS 500)	5.6 x 9.2 x 0.4 (14.2 x 23.4 x 1)
	<b>4000332</b> Sliding canister for 200 mL bag (Gambro DF-200 or equivalent)	6.4 x 7.8 x 0.8 (16.3 x 19.1 x 2)
	<b>4000356</b> Swing arm canister for 200 mL bag (Gambro DF-200 or CryoMACS 250)	6.4 x 7.8 x 0.8 (16.3 x 19.1 x 2)
	<b>4000333</b> Sliding canister for 700 mL bag (Gambro DF-700 or CryoMACS 500 and 750 and OriGen CS 750 and CS 1000)	6.5 x 11.8 x 0.8 (16.5 x 30 x 2)
	<b>4000357</b> Swing arm canister for 700 mL bag (Gambro DF-700 or CryoMACS 500 and 750 and OriGen CS 750 and CS 1000)	6.5 x 12.1 x 0.8 (16.5 x 30.7 x 2)




### Mini square racks for liquid phase storage

Cat. No. and description		Dimensions H x W x D in. (cm)	Storage	CryoPlus 1 system	CryoPlus 2 system	CryoPlus 3 system	CryoPlus 4 system
	4000260	26.2 x 3.3 x 3.4 (66.5 x 8.4 x 8.6)	Racks per vessel	11	30	55	92
	Rack includes cardboard boxes and cardboard dividers for 1.5/2 mL vials		Boxes per rack	12	12	12	12
	Vials per box		25	25	25	25	
	Vials per vessel		3,300	9,000	16,500	27,600	

### Mini boxes for liquid phase storage

	Cat. No.	Description	
	Boxes	189470	Mini cardboard box - 2.75 in. x 2.75 in. x 2 in. for use with 25-cell divider for use with 1.5/2 mL vials
	Cell dividers	189387	25-cell divider for use in the 189470 box



**Racks for Thermo Scientific™ Matrix™ 2D tubes**

Cat. No. and description		Dimensions H x W x D in. (cm)	Storage	CryoPlus 1 system	CryoPlus 2 system	CryoPlus 3 system	CryoPlus 4 system
<b>Liquid phase</b>							
	<b>1950836</b>	27 x 5.313 x 3.65 (68.6 x 13.5 x 9.3)	Racks per vessel	7	16	30	51
	Rack for Matrix 2D plates measuring 1.082 in. to 1.394 in. tall		Plates per rack	18	18	18	18
	Plates per vessel		126	288	900	918	
	<b>1950839</b>	27.15 x 5.313 x 3.65 (69 x 13.5 x 9.3)	Racks Per vessel	7	16	30	51
	Rack for Matrix 2D plates measuring 1.55 in. to 2.15 in. tall		Plates per rack	12	12	12	12
	Plates per vessel		84	192	360	612	
	<b>1950840</b>	27.5 x 5.313 x 3.65 (69.9 x 13.5 x 9.3)	Racks per vessel	7	16	30	51
	Rack for Matrix 2D plates measuring 2.3 in. tall		Plates per rack	11	11	11	11
	Plates per vessel		77	176	330	561	
	<b>1950841</b>	27.5 x 5.313 x 3.65 (69.9 x 13.5 x 9.3)	Racks per vessel	7	16	30	51
	Rack for Matrix 2D plates measuring 3.78 in. tall		Plates per rack	7	7	7	7
	Plates per vessel		49	112	210	357	
<b>Vapor phase</b>							
	<b>1950832</b>	22.5 x 5.313 x 3.65 (57.2 x 13.5 x 9.3)	Racks per vessel	7	16	30	51
	Rack for Matrix 2D plates measuring 1.082 in. to 1.394 in. tall		Plates per rack	15	15	15	15
	Plates per vessel		105	240	450	765	
	<b>1950833</b>	22.6 x 5.313 x 3.65 (57.4 x 13.5 x 9.3)	Racks per vessel	7	16	30	51
	Rack for Matrix 2D plates measuring 1.55 in. to 2.15 in. tall		Plates per rack	10	10	10	10
	Plates per vessel		70	160	300	510	
	<b>1950834</b>	22.5 x 5.313 x 3.65 (57.2 x 13.5 x 9.3)	Racks per vessel	7	16	30	51
	Rack for Matrix 2D plates measuring 2.3 in. tall		Plates per rack	9	9	9	9
	Plates per vessel		63	144	270	459	
	<b>1950835</b>	19.65 x 5.313 x 3.65 (49.9 x 13.5 x 9.3)	Racks per vessel	7	16	30	51
	Rack for Matrix 2D plates measuring 3.78 in. tall		Plates per rack	5	5	5	5
	Plates per vessel		35	80	150	255	


**Platform risers required for racks in vapor phase storage**

	Description of riser and racks for intended use	CryoPlus 1 system	CryoPlus 2 system	CryoPlus 3 system	CryoPlus 4 system
	5 in. riser for use with 1950832, 1950833, 1950834, and 1950835 racks	4000061	4000061	4000064	4000067

### Boxes and dividers

Cat. No.	Description
<b>Boxes</b>	
4000014	Cardboard box with 100-cell divider for use with 1.5/2 mL vials (does have drain hole)
4000015	Cardboard box with 100-cell divider for use with 4/5 mL vials (does have drain hole)
820010	Polycarbonate box with 100-cell divider for use with 1.5/2 mL vials with internal threaded caps (does have drain hole)
820011	Polycarbonate box with 81-cell divider for use with 1.5/2 mL vials (does have drain hole)
4000238	Stainless steel box for 1.5/2 mL vials (does not include 100-cell divider)
4000239	Stainless steel box for 4 mL vials (does not include 100-cell divider)
<b>Cell dividers</b>	
4000013	100-cell divider for use in the 4000014, 4000015, 4000238, and 4000239 boxes

### Other accessories for CryoPlus systems

Cat. No. and description	CryoPlus 1 system	CryoPlus 2 system	CryoPlus 3 system	CryoPlus 4 system	
<b>Recorder kit</b>					
	Factory installed 120 volt, 50/60 Hz with 6 in. diameter circular chart, 7-day temperature recorder, pen and chart drive (Specify top, left and right side mount; right side is default if not specified)		201175		
	Factory installed 220 volt, 50/60 Hz with 6 in. diameter circular chart, 7-day temperature recorder, pen and chart drive (Specify top, left and right side mount; right side is default if not specified)		201280		
	Replacement felt tip pen		245231		
	Replacement circular chart paper		197064		
<b>Hot gas bypass valve assembly</b>					
	Factory installed 120 volt, 50/60 Hz		195021		
	Field installed (QI) 120 volt, 50/60 Hz		195700		
	Factory installed 220 volt, 50/60 Hz		195023		
	Field installed (QI) 220 volt, 50/60 Hz		195701		
<b>Thermal temperature sleeve (27" tall for improved temperature uniformity)</b>					
	Factory installed		Standard, factory-installed option		
		Customer installed/retrofit	1950706	1950707	1950708
<b>Thermal data printed</b>					
	120V, 60 Hz Printer with cable assembly and RS-232 connector kit		4000565		
	220V, 50 Hz Printer with cable assembly and RS-232 connector kit		4000665		
	Thermal printer paper (5 rolls/package)		4000566		
<b>RJ-11 Remote Alarm Cable Kit</b>					
	Used to connect remote alarm to equipment		400142		

### Quality services available (must be ordered with equipment)

260049	Certification of compliance (unit specific)
260045	Temperature mapping
260043	Certificate of conformance