

thermo scientific



TSB Series blood bank refrigerators

Enhanced functionality for secure
blood storage

ThermoFisher
SCIENTIFIC

Security and quality for environmentally sustainable storage of blood and blood components

Thermo Scientific™ TSB Series blood bank refrigerators are designed to meet the strict requirements for storage of blood and blood-related products in full compliance with 93/42/CE—class IIa regulation.



Security

- Digital touchscreen control
- Onboard datalogging tracking temperatures and alarms
- Electronic lock with password protection
- Power failure alarm with battery backup
- Ballasted probe to simulate product temperature




Quality

- 3 year warranty*
- 93/42/CE—class IIa regulation
- White painted interior and precoated steel exterior
- Robust 75 mm thick insulation



Sustainability

- F-gas compliant hydrocarbon green refrigerants
- Low energy consumption with real-time energy monitoring
- Low heat output



*Terms and conditions apply. For more details visit thermofisher.com/refrigerator-services

TSB Series blood bank refrigerators offer:

- A wide range of sizes from 140 L undercounter to 1,500 L double door
- Touchscreen digital controller with onboard datalogging
- Whisper-quiet operation of <50 dBA
- Triple-pane glass doors
- White painted interior and exterior
- Temperature, door ajar, and power failure (with battery backup) alarms
- Low energy consumption with real-time energy monitoring
- Low heat output
- Rear access port
- Electronic lock with password protection

Ordering information

Capacity (L) (450mL Blood Bags)	Doors	Drawers	Dimensions (W x D x H)	Voltage	Plug type	Packaging	Cat. No.
140 (63)	1 glass	3	520 x 635 x 1,000 mm	230 V/50 Hz	EU	Box	TSB140REGBV
						Crate	TSB140REGCV
					UK	Box	TSB140REGBW
						Crate	TSB140REGCW
250 (150)		EU	5		600 x 655 x 1,635 mm	Box	TSB250REGBV
						Crate	TSB250REGCV
		UK	Box		TSB250REGBW		
			Crate		TSB250REGCW		
400 (180)	EU	6	600 x 655 x 1,955 mm	Box	TSB400REGBV		
				Crate	TSB400REGCV		
	UK			Box	TSB400REGBW		
				Crate	TSB400REGCW		
700 (336)	EU	7	750 x 840 x 2,100 mm	Box	TSB700REGBV		
				Crate	TSB700REGCV		
	UK			Box	TSB700REGBW		
				Crate	TSB700REGCW		
1,500 (672)	EU	14	1440 x 825 x 2,070 mm	Box	TSB1500REGBV		
				Crate	TSB1500REGCV		
	UK			Box	TSB1500REGBW		
				Crate	TSB1500REGCW		

Find out more at thermofisher.com/cold